

intellijoint **KNEE**[®] FAQ

Surgeons

Which implant vendor am I required to use with Intellijoint KNEE[®]?

You can use any implant vendor with Intellijoint KNEE. The system is an open platform.

I like my cutting guides; do I have to use Intellijoint specific guides?

You can use your preferred cutting guides with Intellijoint KNEE, but we do provide both an adjustable cutting guide and a mono-block cutting guide as options.

What information does Intellijoint KNEE provide me clinically?

Intellijoint KNEE accurately provides varus / valgus angles for both the distal femoral and the proximal tibial cuts, the flexion angle in the distal femoral cut and the posterior slope of the proximal tibial cut. The system also provides rotational alignment and resection depths for both cuts.

Do I need a CT scan or additional imaging for the Intellijoint KNEE system?

No. There is no imaging needed. The system is registered to the patient through anatomical landmarks including the hip center of rotation, removing the need for an intramedullary (IM) rod.

Are the pins placed outside the surgical incision?

No. Intellijoint KNEE is a pinless system and does not require pins placed outside of your regular incision.

Can I trial the system in the operative room?

Yes. Please contact your local Intellijoint Surgical[®] representative to schedule a clinical evaluation.

How much does the system cost?

Intellijoint KNEE is sold based on a fee-per-use model. The cost includes the Intellijoint KNEE software, disposables, and 24/7 technical support. Please contact your local representative for pricing.

Who is my Intellijoint representative?

To connect you with a representative in your area, please email sales@intellijointsurgical.com or call 1-888-232-2634.

intellijoint KNEE[®] FAQ

Hospital Administrators

A surgeon is interested in Intellijoint KNEE. What is the cost of the system?

Intellijoint Surgical[®] provides Intellijoint KNEE primarily through a fee-per-use model. Capital and operating leases are also available. Please contact your local representative for pricing.

Are there disposables required for the system?

Yes. With a fee-per-use model, the cost includes the technology and the disposables. Intellijoint KNEE is competitively priced.

How much staff training is required?

Minimal staff training is required as Intellijoint HIP is a surgeon-controlled system. Intellijoint will in-service your OR staff prior to an evaluation.

What does Intellijoint KNEE provide the surgeon with?

Intellijoint KNEE provides surgeons with quantitative bone cut and rotational alignment angles, and resection depth measurements to accurately restore the patient's coronal and sagittal alignment using the surgeon's implant of choice. It removes the need for an intramedullary (IM) rod and the use of mechanical alignment jigs.

Do I have to purchase Intellijoint Surgical's implants to use the system?

No, surgeons can continue to use their implants of choice. Intellijoint Surgical does not sell implants.

Do you have data to support improved outcomes?

Yes. Please visit www.intellijoinsurgical.com or contact your local Intellijoint Surgical representative to receive the most up-to-date clinical studies.

Who is my Intellijoint representative?

To connect with a representative in your area, please email sales@intellijoinsurgical.com or call 1-888-232-2634.